

Technische Übersicht

Reziprok 30° Rechtslauf und 150° Linkslauf	Spitzengröße	Max. Durchmesser	16mm	13mm	9mm	6mm	3mm	1mm	Drehmoment g/cm Ncm
Klein 	0.20	1.20	1.20	0.87	0.72	0.54	0.37	0.25	408 – 530 g/cm 4,0 – 5,2 N/cm
Normal 	0.25	1.20	1.20	0.98	0.82	0.61	0.43	0.30	408 – 530 g/cm 4,0 – 5,2 N/cm
Mittel 	0.35	1.20	1.20	0.95	0.82	0.63	0.49	0.40	408 – 530 g/cm 4,0 – 5,2 N/cm
Groß 	0.45	1.20	1.20	0.92	0.82	0.66	0.56	0.50	408 – 530 g/cm 4,0 – 5,2 N/cm