




Technische Übersicht

	Spitzengröße	Max. Durchmesser	16mm	13mm	9mm	6mm	3mm	1mm	Drehmoment (Ncm)	Drehzahl (U/Min)
SX 	0.19	1.20	-	1.20	1.15	0.78	0.40	0.24	2,5 - 3,0	300 - 400
S1 	0.17	1.20	1.20	1.02	0.68	0.49	0.30	0.21	2,5 - 3,0	300 - 400
S2 	0.20	1.20	1.20	1.01	0.77	0.57	0.37	0.25	2,5 - 3,0	300 - 400
F1 	0.20	1.20	1.14	0.98	0.77	0.59	0.42	0.28	2,5 - 3,0	300 - 400
F2 	0.25	1.20	1.15	1.01	0.81	0.65	0.49	0.35	2,5 - 3,0	300 - 400
F3 	0.30	1.20	1.18	1.06	0.89	0.73	0.57	0.41	2,5 - 3,0	300 - 400
F4 	0.40	1.20	1.20	1.06	0.87	0.72	0.57	0.46	2,5 - 3,0	300 - 400
F5 	0.50	1.20	1.20	1.06	0.87	0.75	0.63	0.55	2,5 - 3,0	300 - 400